

Pig CSF3/Granulocyte colony-stimulating factor ELISA Kit



Catalog No: #EK19194

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CSF3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF
Accession No.	O02837
Uniprot	O02837
GeneID	396839;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CSF3

Note: This product is for in vitro research use only