

# Dog CSF2/Granulocyte-macrophage colony-stimulating factor ELISA Kit

Catalog No: #EK19195

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CSF2/Granulocyte-macrophage colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Granulocyte-macrophage colony-stimulating factor
Target Name	Granulocyte-macrophage colony-stimulating factor
Other Names	GM-CSF, Colony-stimulating factor, CSF
Accession No.	P48749
Uniprot	P48749
GeneID	403923;

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CSF2

Note: This product is for in vitro research use only