

Mouse Csf2/Granulocyte-macrophage colony-stimulating factor ELISA Kit

Catalog No: #EK19196

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Csf2/Granulocyte-macrophage colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Granulocyte-macrophage colony-stimulating factor
Target Name	Granulocyte-macrophage colony-stimulating factor
Other Names	GM-CSF, Colony-stimulating factor, CSF, Csfgm
Accession No.	P01587
Uniprot	P01587
GeneID	12981;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100pg/mL Sensitivity: 0.66pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csf2

Note: This product is for in vitro research use only