

Bovine CCL2/C-C motif chemokine 2 ELISA Kit

Catalog No: #EK19197

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CCL2/C-C motif chemokine 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 2
Target Name	C-C motif chemokine 2
Other Names	Acidic seminal fluid protein, Monocyte chemotactic protein 1A, MCP-1, Small-inducible cytokine A2, SCYA2
Accession No.	P28291
Uniprot	P28291
GeneID	281043;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 13.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL2

Note: This product is for in vitro research use only