

Dog CCL2/C-C motif chemokine 2 ELISA Kit

Catalog No: #EK19198

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CCL2/C-C motif chemokine 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog C-C motif chemokine 2
Target Name	C-C motif chemokine 2
Other Names	Monocyte chemoattractant protein 1, Monocyte chemotactic protein 1, MCP-1, Small-inducible cytokine A2, MCP1, SCYA2
Accession No.	P52203
Uniprot	P52203
GeneID	403981;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL2

Note: This product is for in vitro research use only