

Bovine CCL8/C-C motif chemokine 8 ELISA Kit

Catalog No: #EK19199

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CCL8/C-C motif chemokine 8 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 8
Target Name	C-C motif chemokine 8
Other Names	Monocyte chemoattractant protein 2, Monocyte chemotactic protein 2, MCP-2, Small-inducible cytokine A8, MCP2, SCYA8
Accession No.	Q09141
Uniprot	Q09141
GeneID	281044;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL8

Note: This product is for in vitro research use only