

Dog CCL8/C-C motif chemokine 8 ELISA Kit

Catalog No: #EK19200

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CCL8/C-C motif chemokine 8 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog C-C motif chemokine 8
Target Name	C-C motif chemokine 8
Other Names	Monocyte chemoattractant protein 2, Monocyte chemotactic protein 2, MCP-2, Small-inducible cytokine A8, MCP2
Accession No.	Q68AY9
Uniprot	Q68AY9
GeneID	448792;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL8

Note: This product is for in vitro research use only