

Mouse Ccl7/C-C motif chemokine 7 ELISA Kit

Catalog No: #EK19202

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl7/C-C motif chemokine 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 7
Target Name	C-C motif chemokine 7
Other Names	Intercriine/chemokine MARC, Monocyte chemoattractant protein 3, Monocyte chemotactic protein 3, MCP-3, Protein FIC, Small-inducible cytokine A7, Fic, Mcp3, Scya7
Accession No.	Q03366
Uniprot	Q03366
GeneID	20306;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.84pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccl7

Note: This product is for in vitro research use only