

TAS2R4 Antibody

Catalog No: #44628

Package Size: #44628-1 50ul #44628-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAS2R4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TAS2R4 Antibody detects endogenous levels of total TAS2R4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAS2R4
Target Name	TAS2R4
Other Names	TAS2R4, Candidate taste receptor T2R4, T2R4, Taste receptor type 2 member 4
Accession No.	Swiss-Prot#: Q9NYW5NCBI Gene ID:50832
Uniprot	Q9NYW5
GeneID	50832;
Calculated MW	34 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only