

# Mouse Csf1/Macrophage colony-stimulating factor 1 ELISA Kit

Catalog No: #EK19204

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Csf1/Macrophage colony-stimulating factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Macrophage colony-stimulating factor 1
Target Name	Macrophage colony-stimulating factor 1
Other Names	CSF-1,Csfm
Accession No.	P07141
Uniprot	P07141
GeneID	12977;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.09ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Csf1

Note: This product is for in vitro research use only