

Human CCL22/C-C motif chemokine 22 ELISA Kit

Catalog No: #EK19205

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCL22/C-C motif chemokine 22 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 22
Target Name	C-C motif chemokine 22
Other Names	CC chemokine STCP-1,MDC(1-69),Macrophage-derived chemokine,Small-inducible cytokine A22,Stimulated T-cell chemotactic protein 1,MDC,SCYA22,A-152E5.1
Accession No.	O00626
Uniprot	O00626
GeneID	6367;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CCL22

Note: This product is for in vitro research use only