

Mouse Ccl22/C-C motif chemokine 22 ELISA Kit

Catalog No: #EK19206

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Ccl22/C-C motif chemokine 22 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse C-C motif chemokine 22 |
| Target Name | C-C motif chemokine 22 |
| Other Names | Activated B and dendritic cell-derived,CC chemokine ABCD-1,Small-inducible cytokine A22,Abcd1,Scya22 |
| Accession No. | O88430 |
| Uniprot | O88430 |
| GeneID | 20299; |

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.87pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccl22

Note: This product is for in vitro research use only