

Bovine CCL3/C-C motif chemokine 3 ELISA Kit

Catalog No: #EK19207

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CCL3/C-C motif chemokine 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 3
Target Name	C-C motif chemokine 3
Other Names	Macrophage inflammatory protein 1-alpha,MIP-1-alpha,Small-inducible cytokine A3,SCYA3
Accession No.	Q8SQA6
Uniprot	Q8SQA6
GeneID	282170;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCL3

Note: This product is for in vitro research use only