

Bovine CCL4/C-C motif chemokine 4 ELISA Kit

Catalog No: #EK19209

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|------------------------------------------------------|
| Product Name | Bovine CCL4/C-C motif chemokine 4 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Bovine |
| Specificity | Natural and recombinant Bovine C-C motif chemokine 4 |
| Target Name | C-C motif chemokine 4 |
| Other Names | Small-inducible cytokine A4, SCY4A |
| Accession No. | Q17QA1 |
| Uniprot | Q17QA1 |
| GeneID | 414347; |

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL4

Note: This product is for in vitro research use only