

Chicken CCL4/C-C motif chemokine 4 homolog ELISA Kit



Catalog No: #EK19210

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken CCL4/C-C motif chemokine 4 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken C-C motif chemokine 4 homolog
Target Name	C-C motif chemokine 4 homolog
Other Names	Macrophage inflammatory protein 1-beta homolog,MIP-1-beta,Small-inducible cytokine A4 homolog,SCYA4
Accession No.	Q90826
Uniprot	Q90826
GeneID	395551;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCL4

Note: This product is for in vitro research use only