

Dog CCL4/C-C motif chemokine 4 ELISA Kit

Catalog No: #EK19211

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog CCL4/C-C motif chemokine 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog C-C motif chemokine 4
Target Name	C-C motif chemokine 4
Other Names	Small-inducible cytokine A4, SCYA4
Accession No.	Q68AZ0
Uniprot	Q68AZ0
GeneID	448786;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL4

Note: This product is for in vitro research use only