

## Pig CCL4/C-C motif chemokine 4 ELISA Kit

Catalog No: #EK19212

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig CCL4/C-C motif chemokine 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig C-C motif chemokine 4
Target Name	C-C motif chemokine 4
Other Names	Macrophage inflammatory protein 1-beta,MIP-1-beta,Small-inducible cytokine A4,MIP1B,SCY4A
Accession No.	Q711P4
Uniprot	Q711P4
GeneID	396668;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCL4

Note: This product is for in vitro research use only