

TAS2R9 Antibody

Catalog No: #44629

Package Size: #44629-1 50ul #44629-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAS2R9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	TAS2R9 Antibody detects endogenous levels of total TAS2R9
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAS2R9
Target Name	TAS2R9
Other Names	TAS2R9, T2R9, Taste receptor type 2 member 9, TRB6
Accession No.	Swiss-Prot#: Q9NYW1NCBI Gene ID:50835
Uniprot	Q9NYW1
GeneID	50835;
Calculated MW	36 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only