

Rat Ccl4/C-C motif chemokine 4 ELISA Kit

Catalog No: #EK19213

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ccl4/C-C motif chemokine 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat C-C motif chemokine 4
Target Name	C-C motif chemokine 4
Other Names	Macrophage inflammatory protein 1-beta,MIP-1-beta,Small-inducible cytokine A4,Mip1b,Scya4
Accession No.	P50230
Uniprot	P50230
GeneID	116637;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccl4

Note: This product is for in vitro research use only