

## Bovine CCL20/C-C motif chemokine 20 ELISA Kit

Catalog No: #EK19214

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CCL20/C-C motif chemokine 20 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 20
Target Name	C-C motif chemokine 20
Other Names	Macrophage inflammatory protein 3 alpha,MIP-3-alpha,Small-inducible cytokine A20,MIP3A,SCYA20
Accession No.	Q8SQB1
Uniprot	Q8SQB1
GeneID	281666;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCL20

Note: This product is for in vitro research use only