

Mouse Ccl20/C-C motif chemokine 20 ELISA Kit

Catalog No: #EK19215

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl20/C-C motif chemokine 20 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 20
Target Name	C-C motif chemokine 20
Other Names	Beta-chemokine exodus-1, CC chemokine LARC, CC chemokine ST38, Liver and activation-regulated chemokine, Macrophage inflammatory protein 3 alpha, MIP-3-alpha, Small-inducible cytokine A20, Larc, Scya20
Accession No.	O89093
Uniprot	O89093
GeneID	20297;

Application Details

Assay Type: Sandwich Detect Range: 3.9-250pg/mL Sensitivity: 1.97pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccl20

Note: This product is for in vitro research use only