

## Bovine CCL5/C-C motif chemokine 5 ELISA Kit

Catalog No: #EK19218

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Bovine CCL5/C-C motif chemokine 5 ELISA Kit                                   |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Bovine  |
| Specificity        | Natural and recombinant Bovine C-C motif chemokine 5                          |
| Target Name        | C-C motif chemokine 5   |
| Other Names        | SIS-delta, Small-inducible cytokine A5, T-cell-specific protein RANTES, SCYA5 |
| Accession No.      | O97919  |
| Uniprot            | O97919  |
| GeneID             | 327712;   |

## Application Details

Assay Type: Sandwich Detect Range: 39-2500pg/mL Sensitivity: 19.5pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCL5

Note: This product is for in vitro research use only