

Human CCL5/C-C motif chemokine 5 ELISA Kit

Catalog No: #EK19220

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCL5/C-C motif chemokine 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 5
Target Name	C-C motif chemokine 5
Other Names	EoCP, Eosinophil chemotactic cytokine, SIS-delta, Small-inducible cytokine A5, T cell-specific protein P228, TCP228, T-cell-specific protein RANTES, D17S136E, SCYA5
Accession No.	P13501
Uniprot	P13501
GeneID	6352;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL5

Note: This product is for in vitro research use only