

# Mouse Col28a1/Collagen alpha-1(XXVIII) chain ELISA Kit



Catalog No: #EK19223

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Col28a1/Collagen alpha-1(XXVIII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(XXVIII) chain
Target Name	Collagen alpha-1(XXVIII) chain
Other Names	Col28
Accession No.	Q2UY11
Uniprot	Q2UY11
GeneID	213945;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.395ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Col28a1

Note: This product is for in vitro research use only