

Chicken CD40LG/CD40 ligand ELISA Kit

Catalog No: #EK19224

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken CD40LG/CD40 ligand ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken CD40 ligand
Target Name	CD40 ligand
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD40L, TNFSF5
Accession No.	Q9I8D8
Uniprot	Q9I8D8
GeneID	395485;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/ml Sensitivity: 0.088ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CD40LG

Note: This product is for in vitro research use only