

Mouse Cd40lg/CD40 ligand ELISA Kit

Catalog No: #EK19226

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd40lg/CD40 ligand ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CD40 ligand
Target Name	CD40 ligand
Other Names	CD40-L, T-cell antigen Gp39, TNF-related activation protein, TRAP, Tumor necrosis factor ligand superfamily member 5, CD154, Cd40l, Tnfsf5
Accession No.	P27548
Uniprot	P27548
GeneID	21947;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd40lg

Note: This product is for in vitro research use only