

Pig CD40LG/CD40 ligand ELISA Kit

Catalog No: #EK19227

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CD40LG/CD40 ligand ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig CD40 ligand
Target Name	CD40 ligand
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD154, CD40L, TNFSF5
Accession No.	Q95MQ5
Uniprot	Q95MQ5
GeneID	397231;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.092ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CD40LG

Note: This product is for in vitro research use only