

Rat Cd40lg/CD40 ligand ELISA Kit

Catalog No: #EK19228

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cd40lg/CD40 ligand ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CD40 ligand
Target Name	CD40 ligand
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD154, Cd40l, Tnfsf5
Accession No.	Q9Z2V2
Uniprot	Q9Z2V2
GeneID	84349;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.082ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd40lg

Note: This product is for in vitro research use only