

Dog AFP/Alpha-fetoprotein ELISA Kit

Catalog No: #EK19229

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dog AFP/Alpha-fetoprotein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Alpha-fetoprotein
Target Name	Alpha-fetoprotein
Other Names	Alpha-1-fetoprotein,Alpha-fetoglobulin
Accession No.	Q8MJU5
Uniprot	Q8MJU5
GeneID	403551;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.15ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AFP

Note: This product is for in vitro research use only