

Bovine BDH1/D-beta-hydroxybutyrate dehydrogenase, mitochondrial ELISA Kit



Catalog No: #EK19231

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine BDH1/D-beta-hydroxybutyrate dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Target Name	D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Other Names	3-hydroxybutyrate dehydrogenase,BDH,BDH
Accession No.	Q02337
Uniprot	Q02337
GeneID	534090;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BDH1

Note: This product is for in vitro research use only