

Mouse Col4a1/Collagen alpha-1(IV) chain ELISA Kit

Catalog No: #EK19232

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Col4a1/Collagen alpha-1(IV) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(IV) chain
Target Name	Collagen alpha-1(IV) chain
Accession No.	P02463
Uniprot	P02463
GeneID	12826;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.82ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Col4a1

Note: This product is for in vitro research use only