

Human CRYAB/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19233

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRYAB/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin, Heat shock protein beta-5, HspB5, Renal carcinoma antigen NY-REN-27, Rosenthal fiber component, CRYA2, HSPB5
Accession No.	P02511
Uniprot	P02511
GeneID	1410;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CRYAB

Note: This product is for in vitro research use only