

## Mouse Cryab/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19234

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cryab/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin,P23,Crya2
Accession No.	P23927
Uniprot	P23927
GeneID	12955;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cryab

Note: This product is for in vitro research use only