

Rat Cryab/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19235

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cryab/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin
Accession No.	P23928
Uniprot	P23928
GeneID	25420;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cryab

Note: This product is for in vitro research use only