

## Mouse Cav1/Caveolin-1 ELISA Kit

Catalog No: #EK19236

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cav1/Caveolin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caveolin-1
Target Name	Caveolin-1
Other Names	Cav
Accession No.	P49817
Uniprot	P49817
GeneID	12389;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cav1

Note: This product is for in vitro research use only