

Bovine ADM/ADM ELISA Kit

Catalog No: #EK19238



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine ADM/ADM ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine ADM
Target Name	ADM
Accession No.	O62827
Uniprot	O62827
GeneID	280713;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 6.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADM

Note: This product is for in vitro research use only