

## Rat Adm/ADM ELISA Kit

Catalog No: #EK19239

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                 |
|--------------------|---------------------------------|
| Product Name       | Rat Adm/ADM ELISA Kit           |
| Brief Description  | ELISA Kits                      |
| Applications       | ELISA                           |
| Species Reactivity | Rt                              |
| Specificity        | Natural and recombinant Rat ADM |
| Target Name        | ADM                             |
| Accession No.      | P43145                          |
| Uniprot            | P43145                          |
| GeneID             | 25026;                          |

## Application Details

Assay Type: Sandwich Detect Range: 7.8-500pg/mL Sensitivity: 4.17pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adm

Note: This product is for in vitro research use only