

## Mouse Cdo1/Cysteine dioxygenase type 1 ELISA Kit

Catalog No: #EK19242

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdo1/Cysteine dioxygenase type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cysteine dioxygenase type 1
Target Name	Cysteine dioxygenase type 1
Other Names	Cysteine dioxygenase type I,CDO
Accession No.	P60334
Uniprot	P60334
GeneID	12583;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.159ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cdo1

Note: This product is for in vitro research use only