

TAAR6 Antibody

Catalog No: #44632

Package Size: #44632-1 50ul #44632-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TAAR6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TAAR6 Antibody detects endogenous levels of total TAAR6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAAR6
Target Name	TAAR6
Other Names	TAAR6, TA4, TaR-4, TAR4, TAR6, TaR-6, TRAR4, Trace amine receptor 4, Trace amine receptor 6, RP11-295F4.3
Accession No.	Swiss-Prot#: Q96RI8NCBI Gene ID:319100
Uniprot	Q96RI8
GeneID	319100;
Calculated MW	38 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only