

Rat Cdo1/Cysteine dioxygenase type 1 ELISA Kit

Catalog No: #EK19243

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdo1/Cysteine dioxygenase type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cysteine dioxygenase type 1
Target Name	Cysteine dioxygenase type 1
Other Names	Cysteine dioxygenase type I,CDO
Accession No.	P21816
Uniprot	P21816
GeneID	81718;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.16ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdo1

Note: This product is for in vitro research use only