

Mouse Crhr2/Corticotropin-releasing factor receptor 2 ELISA Kit



Catalog No: #EK19245

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Crhr2/Corticotropin-releasing factor receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Corticotropin-releasing factor receptor 2
Target Name	Corticotropin-releasing factor receptor 2
Other Names	CRF-R-2,CRF-RB,Corticotropin-releasing hormone receptor 2,CRH-R-2,Crf2r
Accession No.	Q60748
Uniprot	Q60748
GeneID	12922;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.161ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Crhr2

Note: This product is for in vitro research use only