

Pig CAT/Catalase ELISA Kit

Catalog No: #EK19246

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CAT/Catalase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Catalase
Target Name	Catalase
Accession No.	O62839
Uniprot	O62839
GeneID	397568;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.09ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CAT

Note: This product is for in vitro research use only