

# Mouse CrIs1/Cardiolipin synthase (CMP-forming) ELISA Kit

Catalog No: #EK19247

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Mouse CrIs1/Cardiolipin synthase (CMP-forming) ELISA Kit         |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Ms   |
| Specificity        | Natural and recombinant Mouse Cardiolipin synthase (CMP-forming) |
| Target Name        | Cardiolipin synthase (CMP-forming)                               |
| Other Names        | CLS  |
| Accession No.      | Q80ZM8   |
| Uniprot            | Q80ZM8   |
| GeneID             | 66586;   |

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.408ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CrIs1

Note: This product is for in vitro research use only