

Rat Crls1/Cardiolipin synthase (CMP-forming) ELISA Kit

Catalog No: #EK19248

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Crls1/Cardiolipin synthase (CMP-forming) ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cardiolipin synthase (CMP-forming)
Target Name	Cardiolipin synthase (CMP-forming)
Other Names	CLS
Accession No.	Q5U2V5
Uniprot	Q5U2V5
GeneID	366196;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Crls1

Note: This product is for in vitro research use only