

Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit

Catalog No: #EK19249

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Acetyl-CoA acetyltransferase, mitochondrial |
| Target Name | Acetyl-CoA acetyltransferase, mitochondrial |
| Other Names | Acetoacetyl-CoA thiolase |
| Accession No. | Q8QZT1 |
| Uniprot | Q8QZT1 |
| GeneID | 110446; |

Application Details

Assay Type: Sandwich Detect Range: 7.8-500pg/mL Sensitivity: 4.18pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Acat1

Note: This product is for in vitro research use only