

Rat Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit



Catalog No: #EK19250

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Acetyl-CoA acetyltransferase, mitochondrial
Target Name	Acetyl-CoA acetyltransferase, mitochondrial
Other Names	Acetoacetyl-CoA thiolase
Accession No.	P17764
Uniprot	P17764
GeneID	25014;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.88pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Acat1

Note: This product is for in vitro research use only