

## Rat Apoc1/Apolipoprotein C-I ELISA Kit

Catalog No: #EK19251

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Apoc1/Apolipoprotein C-I ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Apolipoprotein C-I
Target Name	Apolipoprotein C-I
Other Names	Apo-CI, Apolipoprotein C1, Liver regeneration-related protein LRRG04
Accession No.	P19939
Uniprot	P19939
GeneID	108348160;25292;

## Application Details

Assay Type: Sandwich Detect Range: 7.8-500ng/mL Sensitivity: 3.9ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Apoc1

Note: This product is for in vitro research use only