

## Mouse Atoh1/Protein atonal homolog 1 ELISA Kit

Catalog No: #EK19253

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Atoh1/Protein atonal homolog 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Protein atonal homolog 1
Target Name	Protein atonal homolog 1
Other Names	Helix-loop-helix protein mATH-1,mATH1,Ath1
Accession No.	P48985
Uniprot	P48985
GeneID	11921;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.409ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Atoh1

Note: This product is for in vitro research use only