

Mouse Arsa/Arylsulfatase A ELISA Kit

Catalog No: #EK19254

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Arsa/Arylsulfatase A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Arylsulfatase A
Target Name	Arylsulfatase A
Other Names	ASA, Cerebroside-sulfatase, As2
Accession No.	P50428
Uniprot	P50428
GeneID	11883;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.166ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Arsa

Note: This product is for in vitro research use only