

Mouse Anxa5/Annexin A5 ELISA Kit

Catalog No: #EK19255

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Anxa5/Annexin A5 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Annexin A5 |
| Target Name | Annexin A5 |
| Other Names | Anchorin CII,Annexin V,Annexin-5,Calphobindin I,CBP-I,Endonexin II,Lipocortin V,Placental anticoagulant protein 4,PP4,Placental anticoagulant protein I,PAP-I,Thromboplastin inhibitor,Vascular anticoagulant-alpha,VAC-alpha,Anx5 |
| Accession No. | P48036 |
| Uniprot | P48036 |
| GenelD | 11747; |

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Anxa5

Note: This product is for in vitro research use only